


Sani-Cloth®
CLEANING AND DISINFECTANT

Chlor
WIPE



Sani-Cloth® Chlor is a sporicidal cleaning and disinfectant solution ideal for use in high risk, contaminated environments.

- + **Ready to use:** pre-dosed, low linting disposable cleaning and disinfection wipe, requiring no mixing or measuring.
- + **Proven efficacy:** effective against all spores, including C.difficile and fully virucidal in 60 seconds.
- + **Improved efficiencies:** the product solution is impregnated with sodium hypochlorite and surfactant to clean and disinfect providing you with consistently reliable results and significant time savings delivering approximately 6000 ppm available chlorine.
- + **Optimal protection:** can be used on non-invasive medical devices; non-porous hard surfaces with supporting compatibility data.
- + **Cost effective:** using the same solution for cleaning and disinfecting saves you time and money.



+ PDI
BE THE DIFFERENCE™

Product Efficacy

| SANI-CLOTH® CHLOR | TEST METHOD | ORGANISM | CONTACT TIME* | |
|----------------------------|-------------|----------------------------------|---------------|------------|
| Bacteria | EN 16615 | Pseudomonas aeruginosa | 60 sec (c) | 2 min (d) |
| | EN 16615 | Staphylococcus aureus | 60 sec (c) | 2 min (d) |
| | EN 16615 | Enterococcus hirae | 60 sec (c) | 2 min (d) |
| | EN 13727 | Pseudomonas aeruginosa | 30 sec (c) | |
| | EN 13727 | Staphylococcus aureus | 30 sec (c) | |
| | EN 13727 | Enterococcus hirae | 30 sec (c) | |
| Yeasts | EN 16615 | Candida albicans | 60 sec (c) | 2 min (d) |
| | EN 13624 | Candida albicans | 30 sec (c) | |
| | EN 13624 | Candida auris | 30 sec (c) | |
| Enveloped Virus | EN 14476 | Vaccinia Virus | 30 sec (c) | 60 sec (d) |
| | EN 14476 | Corona,HIV, Hep B & C, SARS, Flu | 30 sec (c) | 60 sec (d) |
| Non-enveloped Virus | EN 14476 | Norovirus | 60 sec (c) | 60 sec (d) |
| | EN 14476 | Adenovirus | 60 sec (c) | |
| | EN 14476 | Poliovirus | 60 sec (c) | |
| Spores | EN 16615 | C.difficile | 60 sec (c) | |
| | EN 13704 | Bacillus subtilis | 5 min (c) | |
| | EN 17126 | C.difficile | 15 min (c) | |
| | EN 17126 | Bacillus subtilis | 15 min (c) | |
| | EN 17126 | Bacillus cereus | 15 min (c) | |

*Clean Conditions (c) Dirty Conditions (d)

Point of Use Accessories



Bracket for canisters

Bracket (colour coded): **A001ORG**

Canister brackets affix securely to the wall, holding the canister in place at its location.



Pole Adapter

Product Code: **A002POL**

Pole brackets make it possible to affix the canister brackets on to portable equipment such as BP stands.



I Am Clean Sticky Notes

Product Code: **A005SNT**

Recyclable sticky notes which can be used to easily and effectively identify equipment or surfaces that have been cleaned.



Educational wall charts and signage

Product evaluations, audits, educational wall charts & signage, that can be personalised for your specific needs.

Range Details

| | REORDER NO | WIPE COUNT | WIPE SIZE | QTY PER CASE | CLASSIFICATION |
|---------------------|------------|------------|---------------|--------------|--------------------------|
| 50ct Large Canister | ECH1400 | 50 | 200mm x 195mm | 6 | Medical Device Class IIa |

Other formulations available:



Find out more about how our Layered Approach to Infection Prevention can help you.

PDI International, Pywell Road, Willowbrook Est. Industrial Estate, Corby - Northants, NN17 5XJ, UK

+44 (0) 8081 697 945 contact@pdi-emea.com wearepdi-intl.com PDI International PDI_Intl

