

# Sani-Cloth® 70 is the multi-purpose solution, ideal for disinfection.

- + **Ready to use:** pre-dosed, low linting disposable alcohol disinfection wipes, requiring no mixing or measuring, with no smear formulation.
- + **Proven efficacy:** effective against Coronavirus and Flu in 15 seconds.
- + **Improved efficiencies:** quick drying, no residue formulation containing 70% Isopropanol alcohol allowing the disinfected area or equipment to be used shortly after application.
- + **Maximum flexibility:** can be used on non-invasive medical devices; non-porous hard surfaces.
- + **Recommended usage:** as part of a two-step process with Sani-Cloth® Detergent to clean prior to disinfection.





# **Product Efficacy**

SANI-CLOTH® 70%	TEST METHOD	ORGANISM	CONTACT TIME*	
Bacteria	EN 16615	Pseudomonas aeruginosa	<b>60</b> sec (c)	<b>2</b> min (d)
	EN 16615	Staphylococcus aureus	<b>60</b> sec (c)	2 min (d)
	EN 16615	Enterococcus hirae	<b>60</b> sec (c)	2 min (d)
	EN 13727	Pseudomonas aeruginosa	<b>30</b> sec (c)	
	EN 13727	Staphylococcus aureus	<b>30</b> sec (c)	
	EN 13727	Enterococcus hirae	<b>30</b> sec (c)	
Yeasts	EN 16615	Candida albicans	<b>60</b> sec (c)	2 min (d)
	EN 13624	Candida albicans	<b>30</b> sec (c)	
Enveloped Virus	EN 14476	Vaccinia Virus	<b>15</b> sec (c)	<b>60</b> sec (d)
	EN 14476	Corona,HIV, Hep B & C, SARS, Flu	<b>15</b> sec (c)	<b>60</b> sec (d)
ТВ	EN 14348	Mycobacterium terrae	<b>60</b> sec (c)	

<sup>\*</sup>Clean Conditions (c) Dirty Conditions (d)

# Point of Use Accessories



**Bracket for canisters** Bracket: A001RED

Canister brackets affix securely to the wall, holding the canister in place at its location.



**Pole Adapter** Product Code: A002POL

Pole brackets make it possible to affix the canister brackets on to portable equipment such as BP stands.



### I Am Clean Sticky Notes Product Code: A005SNT

Recyclable sticky notes which can be used to easily and effectively identify equipment or surfaces that have been cleaned.



### **Educational wall charts** and signage

Product evaluations, audits, educational wall charts & signage, that can be personalised for your specific needs.

# Range Details

	REORDER NO	WIPE COUNT	WIPE SIZE	QTY PER CASE	CLASSIFICATION
200ct Large Canister	E701030	200	200mm x 200mm	6	Medical Device Class IIa

# Other formulations available:







## Find out more about how our Layered Approach to Infection Prevention can help you.

PDI EMEA, Pywell Road, Willowbrook Est. Industrial Estate, Corby - Northants, NN17 5XJ, UK











